

LAMPIRAN



2. Pre-treatment dengan APTES
1. Pre-treatment dengan NaOH

Test Protocol http://192.168.1.2/data/NON APPAREL_TREAT2_00154_2022061...

TEXTTEST

INSTRUMENTS

Static Air Permeability

Basic data

Style: NON APPAREL

Reference: TREAT2

Date: 16.06.2022

Time: 18:36:00

Instrument: FX 3300 LabAir IV

Serial Number: 154

Settings

Test pressure: 200 Pa

Test area: 20 cm²

Norm / Min / Max: -1.00 / -1.00 / -1.00 cm³/cm²/s

Statistical analysis

Average: 8.16 cm³/cm²/s

Minimum: 5.13 cm³/cm²/s

Maximum: 12.1 cm³/cm²/s

CV: 32.4 %

Cpk: 0.0000

Test results

1	6.85 cm ³ /cm ² /s	
2	8.53 cm ³ /cm ² /s	
3	5.39 cm ³ /cm ² /s	
4	12.1 cm ³ /cm ² /s	
5	10.9 cm ³ /cm ² /s	
6	5.13 cm ³ /cm ² /s	

1 of 1
6/16/2022 3:49 PM

PM

○ Pengujian Daya Tembus Udara